

PCR Cycler Check™ Advance / OneStep

For conventional PCR block cyclers

INSTRUCTIONS FOR USE

FOR USE IN RESEARCH AND QUALITY CONTROL

Symbols



Lot No.



Cat. No.



Expiry date



Storage temperature



Number of reactions



Manufacturer

INDICATION

False negative PCR results or unspecific amplifications might be caused by a defective PCR cycler. Such cases are critical but can be identified by assessing the temperature accuracy of the PCR cycler. However, temperature assessment of a PCR cycler needs special and therefore expensive equipment, such as temperature sensors that measure the temperature homogeneity in a cycler block.

The PCR Cycler Check™ kit is specifically designed for verifying conventional PCR cyclers, particularly for installation qualification (IQ), operational qualification (OQ) and performance qualification (PQ) as required by various international norms, such as EN ISO 17025, EN 45001, EN ISO 13485, ISO/TS 20836:2007, GLP, GMP, and others.

TEST PRINCIPLE

The PCR Cycler Check™ kit is based on a temperature-sensitive PCR assay to monitor an upper and lower temperature range in one run. The primer sequences in combination with a regular PCR protocol were designed to be extremely sensitive to fluctuations in temperature and thermal homogeneity, precision of the temperature control and timing.

Amplification will be altered when temperature deviates of more than 2 °C from the set value resulting in unexpected band patterns. The cycler performance is tested with typical PCR settings to reflect most users' applications. As an additional indicator of the accurate temperature control of the cycler, the included pre-adjusted target concentrations are only amplified by highly efficient PCRs.

CONTENT

Each kit contains all reagents required to run the PCR. The expiry date of the unopened package is marked on the package label. The kit components must be stored until use at +2 to +8 °C. Do not freeze or store the Validation Reagent after reconstitution.

| Component | Quantity | |
|--------------------|--|--|
| | Advance Cat. No. 57-2102 | OneStep Cat. No. 57-2103 |
| Validation Reagent | Validation Strips: 6 strips, 8 vials each, lyophilized, pre-loaded | Validation Tubes: 4 vials, for 25 reactions each, lyophilized, red cap |
| Caps | 6 cap strips, domed | n.a. |
| Rehydration Buffer | 1 vial, 1.6 ml, blue cap | 2 vials, 1.6 ml each, blue cap |
| Marker | 1 vial, 50 µl, green cap | 2 vials, 50 µl each, green cap |

The lot-specific quality control certificate (Certificate of Analysis) can be downloaded from our website (www.minerva-biolabs.com / www.minervabiolabs.us).

USER-SUPPLIED CONSUMABLES AND EQUIPMENT

The PCR Cycler Check™ kit contains reagents and consumables to perform the cycler check. Additional consumables and equipment are supplied by the user:

- PCR device for 0.2 ml PCR tubes (relevant only for Cat. No. 57-2102).
- Suitable PCR reaction tubes (relevant only for Cat. No. 57-2103)
- 96-well rack for 0.2 ml PCR tubes (relevant only for Cat. No. 57-2102)
- Microcentrifuge for 8-tubes strips (relevant only for Cat. No. 57-2102) and 2 ml reaction tubes
- Vortex
- Pipettes with corresponding filter tips
- Reagents for agarose gel electrophoresis: DNA gel stain, gel running buffer
- Agarose gel electrophoresis equipment and documentation system

PRECAUTIONS


The PCR Cycler Check™ kit is for in vitro use only. The kit should be used by trained laboratory staff only. The PCR Cycler Check™ kit does not contain hazardous substances. Remnants can be discarded according to local regulations.

ADDITIONAL NOTES

- ⇒ These instructions must be understood to successfully use the PCR Cycler Check™ kit. The reagents supplied should not be mixed with reagents from different batches but used as an integral unit. The reagents of the kit must not be used beyond their shelf life.
- ⇒ Follow the exact protocol. Any deviation may affect the test method and can affect the results.
- ⇒ Additional control samples are not required. The kit already contains all necessary controls.

PROCEDURE – OVERVIEW


1A. Reagents preparation for Advance format (Cat. No. 57-2102)




1B. Reagents preparation for OneStep format (Cat. No. 57-2103)



2. PCR amplification



3. Analysis



- Validation Tubes
(only for Cat. No. 57-2103)
- Rehydration Buffer
- Validation Strips
(only for Cat. No. 57-2102)
- Validation Caps
(only for Cat. No. 57-2102)
- Incubate
- Vortex
- Centrifuge
- Ta Annealing Temperature

PROCEDURE - STEP BY STEP

1A. Reagent preparation for Advance format (Cat. No. 57-2102)

1. Spin down the Validation Strips to collect the lyophilized material at the bottom of the tube and place the strips in a 96-well rack. Spin down the rehydration buffer.
 2. Carefully remove the protective seal from the Validation Strips.
 3. Aliquot 25 μ l Rehydration Buffer into each PCR reaction tube. Close the tubes with the provided cap strips.
 4. Incubate for 5 min at room temperature.
 5. Vortex briefly and spin down for 5 sec. Proceed immediately with the PCR.
-


1B. Reagent preparation for OneStep format (Cat. No. 57-2103)

1. Spin down the Validation Tubes and the Rehydration Buffer.
 2. Add 650 μ l of the Rehydration Buffer (blue cap) to each Validation Tube (red cap).
 3. Incubate for 5 min at room temperature.
 4. Vortex briefly and spin down for 5 sec.
-

5. Note: Proceed immediately to step 6. Do not store or freeze the rehydrated Validation Reagent. We recommend reconstituting only the Validation Tube(s) necessary to carry out the selected number of reactions (e.g. 1 vial per 24 reactions, corresponding to 1 cyclor validation).
 6. Aliquot 25 μ l of the rehydrated Validation Reagent into each PCR tube.
 7. Close the PCR tubes and spin down briefly. Proceed immediately with the PCR.
-

2. Perform the PCR cyclers test

Place the PCR tubes in the cycler. We recommend the following scheme depending on the cycler block format:



Program the cycler as follows:

Step 1 (pre-incubation): 94 °C for 2 min

Step 2 (amplification):

Cycles 35

Denaturation 94 °C for 30 sec

Annealing T_a for 30 sec (Annealing Temperature (T_a) is provided on the Certificate of Analysis (CoA))

Elongation 72 °C for 30 sec

Step 3:

Hold 4 °C to 8 °C

3. Analysis

1. Prepare a 1.5 % agarose gel including DNA stain (approx. 5 mm thick, with a 5 mm comb).

Load 5 μ l of each PCR reaction. Load 5 μ l of the provided marker (i.e. customized DNA ladder) in one or more lanes adjacent to the samples lanes.

2. Note: Loading buffer with dye is already included in the mixes. Thus, additional loading buffer or dye is not required.

3. Perform the gel electrophoresis (e.g. 20 min at 100 V).

4. Visualize the PCR results on a suitable transilluminator.

DATA INTERPRETATION

The cycler passed the test if a single band is visible (Fig. 1). The test run is valid but the cycler does not comply with the expected specifications if either no band or two bands are visible.

If no band is visible in any reaction, the experiment should be repeated to exclude a setup mistake. For the re-test, the annealing temperature (T_a) should be reduced by 3 °C to enhance amplification. If the re-test does not show amplification products and the cycler is already suspected to work out of specification, the device should be sent in for service.


If two bands are visible, either the setup of the test was not correct or the cycler is out of specification and should be sent in for service.

Please note, that all PCR reactions must show a uniform result. If this is not the case, most likely one or even more of the Peltier elements have a malfunction. In this case the experiment should be repeated with an adopted loading scheme.

| Fragment size | Interpretation |
|-------------------|--|
| 144 bp and 210 bp | annealing temperature too low denaturation temperature ok |
| 144 bp | Cycler test passed successfully |
| no bands | annealing temperature too high (s. explanation above) or/and denaturation temperature failure |

Fig. 1: Gel figure showing results obtained at different annealing temperatures

- A: Marker
- B: Temperature too low
- C: Temperature correct
- D: Temperature too high



APPENDIX

Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising from of the use, the results of use, or the inability to use this product.

Trademarks

Venor, Onar and AquaScreen are registered trademarks and PCR Cycler Check, ConviFlex, SwabUp, Mycoplasma Off, Food Control, Meat ID, Vegan Control, ExtractNow, WaterShield, LabClean and PCR Clean are trademarks of Minerva Biolabs GmbH, Germany.

Related Products

qPCR Cycler Validation/Qualification

| | | |
|---------|--------------------|---------------|
| 57-2201 | qPCR Cycler Check™ | 100 reactions |
|---------|--------------------|---------------|

PCR Mix

| | | |
|----------------------|--|----------------------|
| 191-0025/-0100/-0250 | ConviFlex™ DNAmix Mix, PCR Mix with Taq polymerase for conventional and qPCR | 25/100/250 reactions |
|----------------------|--|----------------------|

| | | |
|----------------------|---|----------------------|
| 192-0025/-0100/-0250 | ConviFlex™ RT-Taq Mix, RT-PCR Mix with Taq polymerase and retrotranscriptase for conventional and RT-qPCR | 25/100/250 reactions |
|----------------------|---|----------------------|

SwabUp™ Lab Monitoring Kits

| | | |
|----------------|--------------------------------------|---------------|
| 181-0010/-0050 | Sample collection and DNA extraction | 10/50 samples |
|----------------|--------------------------------------|---------------|

Food and Water Assays

| | | |
|----------|--------------------|--------------|
| 36X-X025 | Food Control™ qPCR | 25 reactions |
|----------|--------------------|--------------|

| | | |
|----------------|----------------|------------------|
| 370-1025/-1100 | Meat ID™ Halal | 25/100 reactions |
|----------------|----------------|------------------|

| | | |
|----------------|----------------|------------------|
| 370-2025/-2100 | Vegan Control™ | 25/100 reactions |
|----------------|----------------|------------------|

| | | |
|---------------------|------------------------------------|----------------------|
| 34-2025/-2100/-2250 | AquaScreen® Legionella pneumophila | 25/100/250 reactions |
|---------------------|------------------------------------|----------------------|

| | | |
|---------------------|--------------------------------|----------------------|
| 33-2025/-2100/-2250 | AquaScreen® Legionella species | 25/100/250 reactions |
|---------------------|--------------------------------|----------------------|

| | | |
|---------------------|------------------------------------|----------------------|
| 34-6025/-6100/-6250 | AquaScreen® Pseudomonas aeruginosa | 25/100/250 reactions |
|---------------------|------------------------------------|----------------------|

| | | |
|---------------------|------------------------------|----------------------|
| 34-7025/-7100/-7250 | AquaScreen® Escherichia coli | 25/100/250 reactions |
|---------------------|------------------------------|----------------------|

Contamination Control Kits for conventional PCR

| | | |
|---------------------------|--|-------------------------|
| 11-1025/-1050/-1100/-1250 | Venor®GeM Classic Mycoplasma Detection Kit | 25/50/100/250 reactions |
|---------------------------|--|-------------------------|

| | | |
|---------------------------|--|------------------------|
| 11-7024/-7048/-7096/-7240 | Venor®GeM Advance Mycoplasma Detection Kit | 24/48/96/240 reactions |
|---------------------------|--|------------------------|

| | | |
|---------------------------|--|-------------------------|
| 11-8025/-8050/-8100/-8250 | Venor®GeM OneStep Mycoplasma Detection Kit | 25/50/100/250 reactions |
|---------------------------|--|-------------------------|

| | | |
|---------------------------|------------------------------|-------------------------|
| 12-1025/-1050/-1100/-1250 | Onar® Bacteria Detection Kit | 25/50/100/250 reactions |
|---------------------------|------------------------------|-------------------------|

Contamination Control Kits for qPCR

| | | |
|---------------------|--|----------------------|
| 11-9025/-9100/-9250 | Venor®GeM qEP Mycoplasma Detection Kit | 25/100/250 reactions |
|---------------------|--|----------------------|

| | | |
|------------------------|---|----------------------|
| 11-91025/-91100/-91250 | Venor®GeM qOneStep Mycoplasma Detection Kit | 25/100/250 reactions |
|------------------------|---|----------------------|

Nucleic Acid Extraction

| | | |
|----------------|--------------------------|-------------------|
| 601-1010/-1050 | ExtractNow™ DNA Mini Kit | 10/50 extractions |
|----------------|--------------------------|-------------------|

| | | |
|----------------|--------------------------------|-------------------|
| 602-1010/-1050 | ExtractNow™ Blood DNA Mini Kit | 10/50 extractions |
|----------------|--------------------------------|-------------------|

| | | |
|----------------|--------------------------|-------------------|
| 603-1010/-1050 | ExtractNow™ RNA Mini Kit | 10/50 extractions |
|----------------|--------------------------|-------------------|

| | | |
|----------------|-------------------------|-------------------|
| 604-1010/-1050 | ExtractNow™ CleanUp Kit | 10/50 extractions |
|----------------|-------------------------|-------------------|

| | | |
|----------------|------------------------------|-------------------|
| 605-1010/-1050 | ExtractNow™ Plasmid Mini Kit | 10/50 extractions |
|----------------|------------------------------|-------------------|

| | | |
|----------------|-------------------------------|-------------------|
| 606-1010/-1050 | ExtractNow™ Virus DNA/RNA Kit | 10/50 extractions |
|----------------|-------------------------------|-------------------|

MB Taq DNA Polymerase

| | | |
|---------------------------|--------------------------------|----------------------|
| 53-0050/-0100/-0200/-0250 | MB Taq DNA Polymerase (5 U/μl) | 50/100/200/250 units |
|---------------------------|--------------------------------|----------------------|

| | | |
|---------------------------|--------------------------------|----------------------|
| 53-1050/-1100/-1200/-1250 | MB Taq DNA Polymerase (1 U/μl) | 50/100/200/250 units |
|---------------------------|--------------------------------|----------------------|

PCR Clean™

| | | |
|---------------|--|-------------------|
| 15-2025/-2200 | DNA Decontamination Reagent, Spray bottle/refill bottles | 250 ml/4 × 500 ml |
|---------------|--|-------------------|

| | | |
|---------|---|----------|
| 15-2001 | DNA Decontamination Reagent, Wipes in a dispenser box | 50 wipes |
|---------|---|----------|

| | | |
|---------|---|--------------|
| 15-2002 | DNA Decontamination Reagent, Wipes in refill bags | 5 × 50 wipes |
|---------|---|--------------|

LabClean™

| | | |
|---------|-------------------------------------|-----|
| 15-4100 | DNA Decontamination Reagent, bottle | 1 l |
|---------|-------------------------------------|-----|

WaterShield™

| | | |
|---------------------|--|-----------------------------|
| 15-3015/-3020/-3050 | Water Disinfection Additive for incubators and water baths, 200x concentrate | 15 x 10 ml/3 x 50 ml/500 ml |
|---------------------|--|-----------------------------|

ZellShield™

| | | |
|---------------|---|-----------------|
| 13-0050/-0150 | Contamination Prevention Reagent 100× concentrate | 50 ml/5 x 50 ml |
|---------------|---|-----------------|

Minerva Biolabs GmbH
Schkopauer Ring 13
D-12681 Berlin, Germany

www.minerva-biolabs.com
Ordering: order@minerva-biolabs.com
Support: support@minerva-biolabs.com

USA & Canada

Minerva Biolabs Inc.
1 Jill Ct., Building 16, Unit 10
Hillsborough, NJ 08844
USA

www.minervabiolabs.us
Ordering: order@minervabiolabs.us
Support: help@minervabiolabs.us

Made in Germany

© 2020 Minerva Biolabs
HB18.07EN